



VANTURN

**VANCARE<sup>®</sup>**

We move body and mind



## VanTurn Components:

**A. Upper Handle** This handle is for the caregiver to maneuver and for a patient to use in raising and standing support.

**B. Middle Handle** This handle is for the patient.

**C. Shin Supports** The shin supports can rotate 360 degrees, as well as slide along the horizontal frame.

**D. Right Side Frame Bolt** This bolt, with its convenient turn knob, allows the top frame to be removed from the base for storage or transport.

**E. Outer Caster** There are four outer Casters on the VanTurn that balance the unit.

**F. Center Caster** The center Caster supports the weight and allows for 360 degree turning radius.


**G. Caster Lock** The locking lever with the bright red label allows the caregiver to lock the center Casters for additional stability while raising the patient or during toileting or other functions.

**H. Left Side Frame Bolt** This bolt with its convenient turn knob allows the top frame to be removed from the base for storage or transport.

**I. Feet Placement Decals** The feet placement decals allow the caregiver to assist the patient with the proper placement of their feet.

**J. Stabilizer Bar** By placing their foot here, the caregiver can stabilize the VanTurn as someone raises into a standing position.

**K. Height Adjustment Knob** By loosening this knob and the corresponding knob on the opposite side, you can raise and lower the shin support frame.



## VANTURN FEATURES AND FUNCTIONS

The VanTurn transfer device is used to lift or transfer patients, that still have some strength, to or from beds, wheelchairs or toilets. It is the ideal device in a rehab unit for activating patients. The VanTurn rolls easily with a heavy user and provides a safe and active transfer. As an accessory, the VanTurn belt is available to help the caregiver assist the patient into a standing position, then if desired, secure them into place on the VanTurn. The VanTurn was the first transfer device of its type on the market and offers a new way of safely assisting patients into a sitting to standing position and in managing shorter distance transfers. The product has been developed in close collaboration between users, nursing, staff and designers. Safe patient handling and movement is a high priority in healthcare facilities. This product helps to create a safe working environment for caregivers while delivering a way for the patient to be able to get to a standing position in a natural and dignified way. It is the ideal device in an acute care setting for post-op transfers as well. Available with a patient specific (disposable) belt option for increased cross contamination protection.

## CAR EXTRACTION

Removing a patient from an automobile can be a very risky task. There is little room for a caregiver to maneuver in the open door area. With a VanTurn, a caregiver can safely assist a patient into a standing position, move them away from the auto, and sit them down into a transport chair. It will work equally as well to put a patient back into an automobile.



## VANCARE, INC

1515 1st STREET  
AURORA, NE 68818  
TEL 800.694.4525  
TEL 402.694.4525  
FAX 402.694.3994

[vancare.com](http://vancare.com)

Scan the barcode above with a smart phone for more product information.