


Vera™ II Components

- A. **Model Name** Vera™ II
- B. **Two Point Arm** The 2 point back belt attachment arm.
- C. **Safety Latch** Prevents the loop of the sling From accidentally coming off.
- D. **Adjustable Knee Pad** Stabilizes and comfortably supports.
- E. **Safety Strap** Keeps legs in position and prevents dementia patients from stepping off the foot tray.
- F. **Foot Tray** Heavy-duty plastic foot tray is molded with non-slip tread.
- G. **Caster** 3" casters are hospital grade.
- H. **Locking Caster** The heavy-duty 4" locking casters are hospital grade.
- I. **Scale Load Cells** Scale load-cell housing, only on optional scale models.
- J. **Control Box** Houses electronics including battery gauge, emergency down button and emergency stop button.
- K. **Emergency Stop Button** Used to shut-off power to the lift. Push in to stop power, quarter-turn to restore.
- L. **Acuator** Used to raise and lower the lift arm.
- M. **Manual Emergency Lowering Device** Used to manually lower the lift in the event of power failure.
- N. **Battery Pack** Can be charged on the lift or in a wall mounted charging station
- O. **Manual Base Spread Handle** Used to open/close legs of the lift. Power leg spread option available.
- P. **Hand Control** Used to raise and lower lift as well as to spread the legs with the power leg spread base motor option.
- Q. **Built-in Scale** Digital Scale magnetized for easy display placement.
- R. **Extended Axis Pin** For use with sitting sling and Vera Hygiene sling.
- S. **Handles** Used to maneuver lift into position.



VANCARE®

We move body and mind

Vera™ II
B 350, B 450, B 600



Scan the barcode above with a smart phone for more product information.

VANCARE, INC
1515 1st STREET
AURORA, NE 68818
TEL 800.694.4525
TEL 402.694.4525
FAX 402.694.3994
vancare.com



Made with Pride in the USA

Vera™ II Series

The Vera™ II series of battery-powered lifts allow nursing staff to raise residents/patients from a seated position to a standing position with minimal effort. In most cases, the resident/patient can be transferred from a seated position by one staff member, requiring minimal arm or shoulder assistance from the resident/patient.

The Vera™ II comes in three weight capacities, 350 lbs., 450 lbs., and 600 lbs. Each model can be ordered with a manual base spread or power base spread. Each model can be ordered with a digital scale.

A Vera™ II lift also is a effective and efficient way to reduce injuries among nursing staff, helping to reduce sick leave, workers compensation claims, staff turnover, and liability insurance premiums.



Here's why Vera™ II is appreciated by staff and residents/patients:

- Heavy-duty frame offers tip-resistant support.
- Adjustable base allows for easy access to wheelchairs and geri chairs.
- High quality casters provide easy maneuverability in tight places.
- Dual purpose lift allows for transferring residents/patients via "Sit-to-Stand" transfers for those who are weight bearing and "Sit-to-Sit" transfers for non-weight bearing transfers.

Vera™ II Product Features

- Magnetized hand-held control secures anywhere on lift
- Auxiliary hand-operated control provides means of raising and lowering the resident/patient. In case of power failure, the resident/patient can be safely lowered with an emergency down button. The B450 and B600 model lifts have a manual down/backup capability.
- Removable, and rechargeable 24 volt DC memory-free battery pack features a low-battery power indicator.
- Complies with all applicable federal and state regulations (i.e. OSHA).
- Slip-resistant footplate can be easily removed for convenient cleaning.
- Knee pad can be adjusted for the comfort of the resident/patient.
- Optional Intellipac for keeping track of use and service.
- Optional built in scale.



Safe Patient Handling Options

The Vera™ II Series lifts are a great choice for sit-to-stand lifting needs. But when it comes to Safe Patient Handling, there are more tools necessary to accomplish a safe working environment. That is why Vancare offers a complete Ergonomic Patient Management Program with many equipment options as well as program choices.

