

# VanGo Components

- A. Model Name VanGo.**
- B. Four Point Floating Hanger Assembly** The four point hanger assembly is standard on Vancare lifts. A two-point and four-point clip style is also available.
- C. Safety Latch** Prevents the loop of the sling from accidentally coming off.
- D. Front fork and Castor** The front fork is solid aluminum, no plastic! The dual 3" castor is hospital grade.
- E. Base Socket** This unique base socket design allows for height adjustability and quick disassembly for transporting the lift in a vehicle.
- F. Rear Fork and Locking Castor** The rear fork is solid aluminum, no plastic! The heavy-duty 4" locking caster is hospital grade.
- G. Leg Spread Base Motor** 24v actuator for power spreading of the base legs.
- H. Height Adjust Handle** Use this handle to adjust the height of the mast in three positions.
- I. Control Box** Houses electronics including battery gauge, emergency down button and emergency stop button.
- J. Emergency Stop Button** Used to shut-off power to the lift. Normally in open (out) position.
- K. Hand Control** Used to raise and lower lift as well as to spread the legs with the leg spread base motor.
- L. Battery Pack** Can be charged on the lift or in a wall-mounted charging station.
- M. Height Adjustable Handles** Four height positions.
- N. Actuator** Used to raise and lower the lift arm.
- O. Manual Emergency Lowering Device** Used to manually lower the lift in the event of power failure.
- P. Mast Joint** Cast aluminum mast joint is more durable than plastic.



# VANCARE®

We move body and mind

VanGo  
450, 600



Scan the barcode above with a smart phone for more product information.

**VANCARE, INC**  
1515 1st STREET  
AURORA, NE 68818  
TEL 800.694.4525  
TEL 402.694.4525  
FAX 402.694.3994  
[vancare.com](http://vancare.com)



**Made with Pride in the USA**

## VanGo Features and Functions

The VanGo is the aluminum lift that makes light work of most transfers! The VanGo is a mobile patient lift that can lift patients or residents from a seated or prone position and can easily lift a patient safely from the floor. The VanGo comes in two sizes 450 lbs., and 600 lbs. From small tight spaces to large bariatric patients, there is a VanGo to fit your needs. The raising and lowering functions, as well as the lift's legs, are battery powered for easy positioning with the ergonomically correct four button hand control. In addition to the options listed below a digital scale is an available accessory. The VanGo was designed to be durable in the US healthcare environment. Instead of plastic components, the VanGo uses cast aluminum components for the leg forks and mast joints which eliminate costly replacement parts. The heavy-duty battery pack has the longest battery life on the market. The control box features a battery gauge indicating the power level. A wall charging station is also available.

Safety features include an emergency stop button, manual and electric emergency lowering and safety latches on the hanger assemblies. The VanGo comes standard with the 4 point hanger assembly for the ultimate safe and comfortable lifting!

## VanGo Accessories

Accessories that are available for the VanGo include a 4 Point Stretcher that is used for supine lifting. 2 Point Hanger Assembly is used with 2 Point slings common with European Manufacturers. Ambulation arms are used for gait training and as leverage to assist the caregiver in maneuvering the lift when it is heavy with a patient. 4 point Clip Style Arm is used with clip style slings. Hanging scale for accurately reading patient weights. Intellipack for keeping track of service and use over time.

For more info go to [vancare.com](http://vancare.com)



