


Anatomy of the G-Series Portable



Lifts up to
450 lbs!

Emergency lowering/
Stop feature

Quick magnetic hand
Control connection

Wired hand
Control

Durable all metal gears

G-Series Portable
weighs
6 lbs

Specifications & Listings

Max. Lifting Speed:

No load: 2 inches/second
150lbs: 1.5 inches/second
450lbs: 1 inch/second

Unit weight:

6 lbs

Safe Working Load (SWL):

450 lbs

Dimensions (see diagram):

A: 11" min, 78" max
B: 5.125"
C: 18.125"
D: 19"
Depth of Lift: 3.5"

Number of lifts per charge (Duty : 10/90):

55 with 150 lbs, (25% of strap at midrange)
30 with 450 lbs, (25% of strap at midrange)
Charging time : 2-4 Hours

Lift Case:

Flame retardant ABS

Batteries:

High Capacity, Long Lasting NiMH
1 x 24V (1.6Ah)

Battery Charger:

100-240 VAC 40 Watt

Hand Control:

Wired

Safety:

Emergency stop : Yes
Emergency lowering device : Independent Electrical
Manual Lowering
Upper limit detection
Lower limit detection
Slack tape sensor
Free Fall brake (Over speed governor)
Low battery and Dead battery alarms
Soft start and stop
Overload protection

Approvals:

Can/CSA-C22.2 No 60601-1:08
CAN/CSA-Z10535.1:15

