

Vanturn Dimensions			
Name	Vanturn Short	Vanturn	Vanturn 450
Part Number	4-1-7400	4-1-7500	4-1-7600
Unit	1 piece	1 piece	1 piece
Baseplate Length	27 in	27 in	31.3 in
Baseplate Width	22.4 in	22.4 in	26.5 in
Baseplate Height	1.6 in	1.6 in	1.8 in
Overall Height	39.2 in	45.3 in	47.6 in
Safe Working Load	330 lbs	330 lbs	452 lbs
Total Weight of lift	36.4 lbs	36.4 lbs	57.3 lbs