

VERA LIFT SPECIFICATIONS

MODEL 350

CONSTRUCTION

- * Smooth, electric worm-gear lift drive (not hydraulic, chain, or belt driven).
- * Hand-held controller (push-button operated, heavy duty coil cord, operable with single finger touch, magnetized so as to secure anywhere on the lifter frame).
- * Back-up safety switches on lift base (up/down control and manual shut-off).
- * 3" dual front casters, 4" single rear locking casters - Standard (low and high options available)
- * Fully-padded, patient contact points.
- * Welded steel (no aluminum); chip, peel, and rust proof hammertone finish; no painted or plastic coatings on steel components.
- * 24 volt, rechargeable "memory free" batteries, gel-cel (not lead acid)
- * Built in battery charger (in the lift base) with overcharging protection; approximately 125-150 transfers per charge; recharge time (15 minute quick charge) and 6-8 hours full charge.

LIFT

- * Fixed mast with pivoting boom.
 - * Expandable base; hand-lever activated to accommodate bed, chair, toilet, or most specialized equipment.
 - * Capacity to clear all under-bed frames.
 - * Height-adjustable foot and knee rest with leg support strap.
 - * Large lifting capacity, up to 350 pounds.
 - * Ability to maneuver through small doorways and bathrooms.
- Lift weight: 138 lbs

BACK BELTS

- * Detach from the lifter when transfer is complete.
- * Accommodate 50 to 350 pound residents (three sizes of belts).
- * Allow toileting without soiling the sling.
- * Use no chains or metal support.
- * Provide safety strap for patient's chest with easy press-to-release buckle.
- * Are machine washable (up to 180° C, quick drying; conveniently stored in pouches).

DIMENSIONS

- | | |
|--------------------------------------|------------------------------------|
| * Overall width with legs closed 25" | * Overall width with legs open 36" |
| * Chassis width 24" | * Chassis height - 4 1/4" |
| * Chassis length - 42 1/2" | * Overall length 43 1/2" |
| * Casters - 3" wheels (4" optional) | * Length of boom 32" |
| * Height to top of tower - 51" | |