



CEILING LIFT ARMOIRE

VANCARE®


We move body and mind



Don't like the looks of a ceiling lift sitting in the corner of your rooms? Bring the look of home into your facility with the Vancare Ceiling Lift Armoire to discretely hide your ceiling lift.

Features:


- Easily & securely attaches to any wall
- Armoire doors have lock & key as standard feature
- Designed to work with straight or xy rail systems
- Easy to clean floor under cabinet
- Edges sealed for moisture barrier & bacteria resistance
- Shelves are height adjustable
- Storage for safe patient handling accessories
- Can be customized to fit the design requirements
- Available in 12 standard finishes
- Saves room over traditional on the floor furniture
- Many styles available to match existing furniture
- Locked Armoire meets the needs of cognitively impaired environments
- Made in the USA



* With Drawers

FEATURES AND FUNCTIONS

The Ceiling Lift Armoire is designed to work with straight rails, as well as with traversing rail systems. It is installed on the wall below the ceiling lift motor. The bi-fold cabinet doors can be locked for securing the motor. The top shelf has a cutout to allow the hanger assembly to be stored at a lower level that will be easily accessible by staff. The bottom two drawers can be used for clothing or to store additional slings and accessories.



* Without Drawers

ON SITE EVALUATION

Every ceiling lift project has it's own criteria to achieve the desired results. It is important to remember that your needs today may not be your needs two years from today. Discuss options with your Vancare representative to be assured that you will get the most value for your dollars spent. Having an on-site evaluation will put you on the right track!



VANCARE, INC
1515 1st STREET
AURORA, NE 68818
TEL 800.694.4525
TEL 402.694.4525
FAX 402.694.3994
vancare.com