

VanTurn Dimensions

Name	Vanturn Short	Vanturn	Vanturn 450
Part Number	7400	7500i	7600
Unit	1 piece	1 piece	1 piece
Baseplate Length	27in	27in	31.3in
Baseplate Width	22.4in	22.4in	26.5in
Baseplate Height	1.6in	1.6in	1.8in
Overall Height	39.2in	45.3in	47.6in
Safe Working Load	330lbs	330lbs	452lbs
Total weight of lift	36.4lbs	36.4lbs	57.3lbs