

LIFT MODELS	C-300	C-450	C-625	C-1000	P-300	P440	P-600
Fixed	•	•	•	•			
Portable					•	•	•
Manual Traverse	•	•	•	•	•	•	•
Power Traverse	•	•	•	•			
Power X-Y System	•	•	•	•			
Lifting weight capacity	300 lbs (136 kg)	450 lbs (204 kg)	625 lbs (283 kg)	1000 lbs (454 kg)	300 lbs (136 kg)	440 lbs (200 kg)	600 lbs (273 kg)
FEATURES							
Pneumatic controls	•	•	•	•	•	•	•
On-board controls					•	•	•
Digital display (LCD)	•	•	•	•			•
Audible low battery	•	•	•	•	•	•	•
Soft start-stop	•	•	•	•	•	•	•
Auto shut-off	•	•	•	•	•	•	•
Changeable carry bar	•	•	•	•			
Return to charge option	•	•	•				
SAFETY FEATURES							
Emergency stop	•	•	•	•	•	•	•
Emergency power lowering	•	•	•	•	•	•	•
Emergency manual raising or lowering	•	•	•	•	•		•
Fail safe over-speed governor	•	•	•	•	•	•	•
Built-in safety limits	•	•	•	•	•	•	•
Slack strap sensor	•	•	•	•	•	•	•
SPECIFICATIONS							
24 volt battery system	•	•	•	•	•	•	•
Up/down speed	1.8"/sec (4.5 cm/sec)	1.8"/sec (4.5 cm/sec)	1.8"/sec (4.5 cm/sec)	.6"/sec (1.5cm/sec)	1.8"/sec (4.5 cm/sec)	1.5"/sec (3.75 cm/sec)	1.8"/sec (4.5 cm/sec)
Traverse speed		4 programmable speeds	4 programmable speeds	4 programmable speeds			
Lift weight	5.7 lbs (2.6 kg)	22 lbs (10 kg)	22 lbs (10 kg)	24.5 lbs (11.1 kg)	7 lbs (3.2 kg)	10.5 lbs (4.8 kg)	24 lbs (10.9 kg)
Lifting height	77" (1.96 m)	80" (2.0 m)	80" (2.0 m)	96" (2.43 m)	77" (1.96 m)	80" (2.0 m)	80" (2.0 m)